



Quality. Service. Value.

[Go to www.tps.com.au](http://www.tps.com.au) for the latest information on the complete TPS product range

TPS has been certified to ISO9001:2000



[Click to go to the Quality Assurance section of the TPS Web Site.](#)

TPS is proud to announce that we have now been certified to the ISO9001:2000 Quality Assurance standard by SGS Systems & Services Certification.

TPS had previously been certified to ISO9001:1994, and the independent auditor made a special mention of the fact that our Quality Management System needed very few changes to meet the new standard.

This was largely due to customer and service focus that has always been a hallmark of TPS.

A copy of our ISO9001:2000 certificate may be downloaded from the [Quality Assurance](#) section of the TPS Web Site.

New Cordwheel cable reel now available



[Click to go to directly the Cordwheel product page.](#)

Part Number 130042

The Cordwheel cable reel is ideal for users of portable instruments which are fitted with long cables.

The Cordwheel is lightweight and easy to use. It can be set up for right or left handed operation.

The Cordwheel can hold 50 or more metres of cable, depending on cable diameter.

The recommended selling price for the Cordwheel is \$39 + GST.

New Flow Cell for laboratory use



[Click to go to directly the Flow Cell product page.](#)

Part Number 112303

Many customers have asked for flow-through Conductivity, pH, ORP and Ion Selective electrodes. In the past, this has meant having a separate version of each and every type of sensor in a flow-through configuration.

Now TPS has introduced a new Flow Cell for laboratory use. It is universally adaptable to all 12mm diameter sensors. Any standard sensor can now be used in a flow-through application.

The recommended selling price for the Flow Cell is \$155 + GST.

Simple, low cost probe stand



[Click to go to directly the Probe Stand product page.](#)

Part Number 121344

TPS has had this Probe Stand in its product range for a number of years. As we have had a number of customers with portable instruments asking for this type of product, we thought it a good idea to mention it in this newsletter.

The stainless steel rod easily unscrews from the base, making the Probe Stand very portable. It weighs only 150 grams and will fit into most carry cases.

At \$28 + GST recommended selling price, it's great value.

TPS Pty Ltd

4 Jamberoo St, Springwood,
Brisbane, Australia, 4127

Phone Australia (07) 32 900
400

International +61 32 900
7 400

Fax Australia (07) 3808 4871

International +61 3808 4871
7

Email tps@tps.com.au

Web www.tps.com.au

Date of Issue : 10-Sep-2003.



[Download a printer friendly PDF version now.](#)

Hint : Right click on the link, then select "Save Target As..." or "Save Link As..." to save the file to your own computer.

If you do not have Adobe Acrobat in your computer, a free Acrobat Reader can be downloaded by clicking on the following link.

