



Quality. Service. Value.

[TPS Home](#)[Contact Us](#)[Products](#)[Distributors](#)[Prices](#)

Flexible Probe Arm



请使用翻译作用于这页底部以翻译页对汉语。

TPS Australia - Quality hand-made instruments since 1968.

Flexible Arm Style Probe Holder



(sensor & beaker not included)

- ◆ Holds 2 x 12mm and 1 x 6.5mm sensor
- ◆ One handed operation
- ◆ Heavy base to stop holder tipping with multiple sensors
- ◆ 1.7 kg shipping weight
- ◆ Order No. 130088

[Full Size Picture](#)



[Download this Specification Sheet in PDF format.](#)

Hint : Right click on a link, then select "Save Target As..." or "Save Link As..." to save the file to your own computer.

If you do not have Adobe Acrobat in your computer, a free Acrobat Reader can be downloaded by clicking on the following link.



TPS reserves the right to change any part of this specification without notice.
Version 1.0
Copyright © 2003 T.P.S. Pty Ltd

[WP-Sim Simulator](#)[Instrument Carry Case](#)[90 Series Case](#)[Flexible Probe Arm](#)

Send mail to tps@tps.com.au with questions or comments about this web site.
Web Author : Michael Schimkat
Copyright © 2002, T.P.S. Pty Ltd
Last modified: April 17, 2003

