

[TPS Home](#)[Contact Us](#)[Products](#)[Distributors](#)[Brochures](#)[Prices](#)

Serial to USB Adaptor



请使用翻译作用于这页底部以翻译页对汉语。

TPS Australia - Quality hand-made instruments since 1968.

XAUSB Serial to USB Adaptor



- ◆ *Full compliance with USB Spec. Rev 1.1*
- ◆ *Win98, Win98 SE, Win ME, Win2000 & Windows XP Compatible*
- ◆ *iMac compatible*
- ◆ *200g shipping weight*
- ◆ *Order No. 130087*

The **XAUSB** USB Serial to USB Adaptor allows TPS instruments with Serial (RS232) interfaces to be connected to a USB port on a Windows or iMac computer. The **XAUSB** is the ideal solution for Laptop or Notebook computers, which often are not supplied without a serial port.

Once the **XAUSB** has been connected and the driver has been correctly installed, the computer will appear to have an extra COM port. Communication software, such as [WinTPS](#), can then connect to your TPS instrument via this new COM port.



[Download this Specification Sheet in PDF format.](#)



[Download the Windows 98, ME, 2000 & XP Drivers for the XAUSB.](#)



[Download the Mac-OS Drivers for the XAUSB](#)

Hint : Right click on a link, then select "Save Target As..." or "Save Link As..." to save the file to your own computer.

If you do not have Adobe Acrobat in your computer, a free Acrobat Reader can be downloaded by clicking on the following link.



*TPS reserves the right to change any part of this specification without notice.
Version 1.0
Copyright © 2003 T.P.S. Pty Ltd*

[WP-Sim Simulator](#)

[Instrument Carry Case](#)

[90 Series Case](#)

[Flexible Probe](#)

[Serial to USB Adaptor](#)

Send mail to tps@tps.com.au with questions or comments about this web site.

Web Author : Michael Schimkat

Copyright © 2002, T.P.S. Pty Ltd

Last modified: June 26, 2003

