



Quality. Service. Value.

- TPS Home
- Contact Us
- Distributors
- Prices

Sensor Connectors



请使用翻译作用于这页底部以翻译页对汉语。

TPS Australia - Quality hand-made instruments since 1968.

Connectors for pH, Redox (ORP), Ion and Reference Sensors




Type : BNC0 Part No : 121402

For all TPS meters, other than the 1852 series and some other older models. Many other manufacturers use the BNC connector.



Type : TPS0 Part No : 121406

For mono pH, Redox, Ion or metal sensors on TPS 1852 & 1952 Series and models 803 and 804. Also for Tacussel and Lyons (Solea) brands.



Type : TPS4 Part No : 121404

For combination pH, Redox or Ion sensors on TPS 1852 & 1952 Series and models 803 and 804. Also for Tacussel and Lyons (Solea) brands. Supplied complete with REF4 reference connector.



Type : COR0 Part No : 121418

For Corning portable meters. May also need REF2 connector.



Type : USA0 Part No : 121410

For mono pH, Redox, Ion or metals sensors on Beckman, Corning, Coleman, Fisher, Horiba, IL, Leeds & Northrup, Orion, Photovolt and Sargeant brands.

Type : USA2 Part No : 121408

For combination pH, Redox, Ion or metals sensors on Beckman, Corning, Coleman,



Fisher, Horiba, IL, Leeds & Northrup, Orion, Photovolt and Sargeant brands. Supplied complete with REF2 connector.



Type : RAS0 **Part No : 121414**

For Radiometer and Crison brands. May also need REF2 connector, below.



Type : REF4 **Part No : 121453**

4mm diameter connector for reference sensor on all TPS meters which are fitted with a separate reference socket. Also used by many other manufacturers.



Type : REF3 **Part No : 121452**

3mm diameter connector for reference sensor on EIL, Pye, Pye-Unicam, Philips and Jones brands.



Type : REF2 **Part No : 121450**

2mm diameter pin connector for reference sensor on Accurex, Beckman, Chemtrix, Corning, Coleman, CHK, Fisher, Horiba, IL, Jenco, Jenway, Kane-May, Leeds & Northrup, Orion, Portec, Photovolt, Sargeant, Starna and Watson-Victor brands.



Type : BELO **Part No : 121424**

TPS connector for model HP80-A. Also for Anax, Activon, EIL, Jenway, Pye, Pye-Unicam, WPA and Philips brands.



Type : BEL3 **Part No : 121422**

For EIL, Pye-Unicam, Pye and Philips brands. Supplied with REF3 connector.



Type : NLJ3 **Part No : 121426**


For N.L. Jones brand meters only.



Type : MET0 **Part No : 121412**

DIN 'E' connector. Also used for Behr, Consort, Decma, DDS, IL, Mettler, Metrohm, Polymetron, Schott-Gerate, Toptronic, and

WTW brands.




Type : RAL4 **Part No : 121416**

For early Radiometer meters. Supplied complete with REF4 reference connector.



Type : HOR0 **Part No : 121420**

For Philips, Elektrofakt and Horiba brands. May also need REF4 reference connector.



Type : PHL4 **Part No : 121428**

For Philips, Elektrofakt and Horiba brands. May also need REF4 reference connector.



Type : ANA0 **Part No : 121432**

For later model Metrohm meters.



[Download this Specification Sheet in PDF format.](#)

Hint : Right click on a link, then select "Save Target As..." or "Save Link As..." to save the file to your own computer.

If you do not have Adobe Acrobat in your computer, a free Acrobat Reader can be downloaded by clicking on the following link.



*TPS reserves the right to change any part of this specification without notice.
Version 1.0
Copyright © 2003 T.P.S. Pty Ltd*

- [pH Sensors](#)
- [Redox Sensors](#)
- [Specific Ion Sensors](#)
- [Conductivity-TDS Sensors](#)
- [Dissolved Oxygen Sensors](#)
- [Temperature Sensors](#)
- [Turbidity Sensors](#)
- [Other Sensors](#)
- [Sensor Connectors](#)

Send mail to tps@tps.com.au with questions or comments about this web site.
Web Author : Michael Schimkat
Copyright © 2002, T.P.S. Pty Ltd
Last modified: January 20, 2003

