



Quality. Service. Value.

TPS Home

Contact Us


Distributors

Prices

WinTPS

TPS Australia - Quality hand-made instruments since 1968.

WinTPS RS232 Communication Software for TPS Instruments



- ◆ *Universal communication software for all TPS instruments with RS232 ports.*
- ◆ *Easy to use with large icons on toolbar.*
- ◆ *Complete help files.*
- ◆ *Compatible with Microsoft™ Windows™ 95, 98, ME, 2000, XP & NT*

[Larger screen shot](#)

WinTPS provides the following functions...

- ◆ Download current data from instrument.
- ◆ Download data stored in instrument's log memory.
- ◆ Download GLP calibration data.
- ◆ "Listen" mode waits data to be sent from the instrument to the computer.
- ◆ Erase instrument's memory remotely.
- ◆ Enable and disable automatic datalogging for 90 series instruments.
- ◆ Start and stop the dissolved oxygen stirrer for 90 and smartCHEM series instruments with Dissolved Oxygen readout.
- ◆ Remotely switch 90 series instruments ON and OFF.

Ordering Information

Part No

WinTPS Communication.....130086

Want more information ?

Download the handbook below for complete details.



[Download this Specification Sheet in PDF format.](#)



[Download the WinTPS handbook.](#)

Hint : Right click on a link, then select "Save Target As..." or "Save Link Ss..." to save the file to your own computer.

If you do not have Adobe Acrobat in your computer, a free Acrobat Reader can be downloaded by clicking on the following link.



TPS reserves the right to change any part of this specification without notice.

Version 1.0

Copyright © 2002 T.P.S. Pty Ltd

Send mail to tps@tps.com.au with questions or comments about this web site.

Web Author : Michael Schimkat

Copyright © 2002, T.P.S. Pty Ltd

Last modified: October 23, 2002